

FLT3 ITD MRD Assay

RUO Kit - Now Available for MiSeq™



BACKGROUND

Stratifying acute myeloid leukemia (AML) disease according to molecular genetic alterations, such as those in the fms related tyrosine kinase 3 (*FLT3*) gene, aids in prognosis¹. The most frequent and clinically significant type of *FLT3* mutation is an internal tandem duplication (ITD) in the juxtamembrane domain². The *FLT3*-ITD mutation occurs in about 25% of newly diagnosed AML patients and is associated with an increased risk of relapse and lower overall survival rate¹. Unlike flow cytometry assays which require fresh sample and are highly subjective, the *FLT3* ITD MRD Assay, a targeted, deep sequencing assay can be used with previously isolated DNA to detect ITD mutations at an allelic sensitivity level of 5×10^{-5} . To further simplify your workflow, our *FLT3* ITD MRD v1.2 software automates data analysis and provides an objective variant call to easily monitor and track efficacy-response assessments and to further guide treatment.

Invivoscribe's *FLT3* ITD MRD Assay includes 24 unique dual-indices enabling the ability to multiplex multiple samples. This kit configuration provides laboratories the flexibility to scale testing for variable AML MRD research needs.


KEY BENEFITS

- » Bring MRD testing in-house for faster, cost-effective results
- » Streamlined workflow reduces errors
- » Flexibility to multiplex samples
- » Scalable for low and high throughput labs
- » Use previously isolated DNA enabling sample batching
- » Dockerized software enables highly portable, flexible and efficient sample analysis
- » Standardized for international collaboration

PRINCIPLE OF THE PROCEDURE

FLT3 ITD or length mutations are caused by duplication and insertion of a portion of the *FLT3* gene that includes the region in and around the juxtamembrane (JM) region. Next-generation sequencing of the PCR products is used to identify DNA sequences specific to previously identified mutations and estimate variant read frequencies (VRF). The software, *FLT3* ITD MRD v1.2, provides an objective variant call in a .tsv output file to automate your AML studies.

WORKFLOW



For Research Use Only (RUO). Not intended for diagnostic purposes.

ORDERING INFORMATION

Catalog #	Product	Quantity
14120019	<i>FLT3</i> ITD MRD Assay (MiSeq™)	96 reactions
14120029	<i>FLT3</i> ITD MRD v1.2 Software (MiSeq™)	1 Dockerized Application

REAGENTS INCLUDED IN THE KIT

Controls	Quantity
<i>FLT3</i> ITD Positive Control	500 µL tube x 2 each
<i>FLT3</i> ITD Negative Control	500 µL tube x 2 each

Master Mixes	Quantity
<i>FLT3</i> ITD Master Mixes	75 µL tube x 24 each

VAF READ DEPTH

FLT3 ITD mutated subjects enrolled in the CHRYSALIS study, who were treated with *FLT3*-inhibitory oral doses of 120mg/day or 200mg/day gilteritinib, had their molecular response assessed from bone marrow aspirates obtained at baseline and at ≥1 additional time point. *FLT3* ITD and total *FLT3* alleles were quantified using the Invivoscribe® *FLT3* ITD MRD assay and used to determine molecular response³. A Cox regression model of overall survival (OS) by Kaplan-Meier estimation was used to evaluate the impact of ITD variant allele frequency (VAF) on overall survival. Molecular response was defined as follows:


Molecular response = ITD VAF (*FLT3* mutant reads: *FLT3* total reads) of $\leq 10^{-2}$ point.

Major molecular response = ITD VAF of $\leq 10^{-3}$

Negative MRD status = ITD VAF of $\leq 10^{-4}$

As shown in Table 1.0 and Figure 1.0, subjects with molecular response had longer overall survival than those without a molecular response³. This demonstrates that when evaluating MRD, achieving a sensitivity $> 10^{-2}$ may be unnecessary when evaluating VAF for OS.

Figure 1.0



Mark J. Levis, et al. (2018) *Blood Adv.* 2(8):825-831.

Table 1.0

Molecular response	Achieved a molecular response		Did not achieve a molecular response		P
	n	Median OS (95% CI), d	n	Median OS (95% CI), d	
ITD VAF $\leq 10^{-2}$	20	417 (246-NA)	60	199 (142-234)	<.001
ITD VAF $\leq 10^{-3}$	18	417 (228-NA)	62	213 (143-264)	.003
ITD VAF $\leq 10^{-4}$ (MRD negative)	13	417 (228-NA)	67	213 (144-264)	.002

Table 1.0 Overall Survival of AML Subjects with *FLT3* ITD

Comparison between subjects achieving a molecular response (*FLT3* ITD VAF $> 10^{-2}$) by the MRD assay and those not achieving a molecular response by the MRD assay. The P values were determined by the log-rank test.

Figure 1.0 Subjects, Overall Survival Stratified by Molecular Response, using the Internationally-Harmonized *FLT3* ITD MRD Assay.

REFERENCES

1. Naval Daver, et al. (2019) *Leukemia.* 33:299-312.
2. Heiko Konig, Mark Levis, (2015) *Expert Opin. Ther. Targets.* 19(1): 37-54
3. Mark J. Levis, et al. (2018) *Blood Adv.* 2(8):825-831.

